**MINI PEDALIER
PEDAL EXERCICER
434017**

Description – Description :

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The mini footswitch is a perfect tool to tone your arms and legs in your home. Compact and portable, this mini footswitch can be easily placed in different rooms or closets in your home, or even under a table or desk to do your exercises while sitting. Lightweight, you can also place it on a table to tone and strengthen your arms. An adjustment knob will allow you to adjust the resistance at several levels, according to your requirements. Sturdy, it is made of sturdy steel.

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4. Head of bolt should be properly positioned into pre-cut slot on leg brace. Using a wrench or pliers, tighten the nuts from the bottom until secure.

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5. Slide the four non-skid tips (G) onto the ends of the leg braces until secure. Using a small amount of liquid soap will assist in applying tips.

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6. Check pedal exerciser for balance by making sure that all four non-skid tips touch the floor surface when exerciser is positioned on the floor surface.

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7. If the exerciser is not balanced, (1) check the bolts for tightness, (2) check the alignment of the leg braces, (3) check the positioning of the non-skid tips.

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8. Insert tension knob (F) and tighten as desired.

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Tighten all bolts periodically. Bolts may require occasional tightening after use.

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The metal components on the tension knob may become hot when in use.

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This product can be used to exercise the legs while sitting in a chair. Place product on the floor at a comfortable distance from the chair – not too far away to cause stress or overstretching of the legs.

The arms can be exercised by placing the product on a table. User sits in a chair at the table and grasps with hands and rotates pedals with hands.

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Please inspect all parts carefully for cracks/breaks prior to use.

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DO NOT STAND ON THIS PRODUCT. IT IS DESIGNED FOR USE WHILE SITTING IN A CHAIR.

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Persons with severe disabilities or persons who are unable to walk without an aide should not attempt to use this product without assistance.

INSTRUCTIONS D'ENTRETIEN - CARE INSTRUCTIONS

Wipe with a clean cloth. If product becomes wet, dry immediately with towel to avoid rust.

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