VERMEIREN

Albatros II

INSTRUCTION MANUAL
MODE D'EMPLOI
GEBRUIKSAANWIJZING
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUKCJA OBSŁUGI

EN

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Contents

| | |
|--|-----------|
| Preface | 2 |
| 1 Product description | 3 |
| 1.1 Intended Use..... | 3 |
| 1.2 Technical specifications | 4 |
| 1.3 Drawing | 6 |
| 1.4 Explanation of the symbols | 7 |
| 1.5 Delivery..... | 7 |
| 2 Use | 7 |
| 2.1 General notes..... | 7 |
| 2.2 Use of the JUMBO system (battery, control box + included battery charger)..... | 10 |
| 2.3 Operating the lift | 11 |
| 2.4 Slings..... | 14 |
| 2.5 Safety rules..... | 14 |
| 2.6 Emergency..... | 15 |
| 3 Installation and adjustment | 16 |
| 3.1 Swapping the battery | 16 |
| 3.2 Assembly or dismantling..... | 16 |
| 4 Maintenance | 23 |
| 5 Ordernumbers | 23 |

Preface

First of all we want to thank you for putting your trust in us by selecting one of our products.

The Vermeiren patient hoists are the result of many years of research and experience. During the development, special attention was given to the ease of use and the serviceability of the patient hoist.

The expected lifetime of your patient hoist is strongly influenced by the care and maintenance of the patient hoist.

This manual will help you get acquainted with the operation of your product.

Following the user instructions and the maintenance instructions are an essential part of the guarantee.

This manual reflects the latest product developments. Vermeiren has the right to introduce changes without the obligation to adapt or replace previously delivered models.

For any further questions, please consult your specialist dealer.

1 Product description

1.1 *Intended Use*

The patient hoist is intended for people with walking difficulties or no walking abilities.

The patient hoist is designed to transport 1 person.

The patient hoist is suited for indoor use.

The user must be lifted in the patient hoist by an attendant.

The different possibilities/positions allow the application of the patient lift in case of inability or impossibility to walk based on:

- paresis
- limb deformations
- cachexia
- atrophy
- and also for aged persons.

When providing for individual requirements:

- body size and weight (max. 150 kg)
- physical and psychological condition
- residential circumstances
- environment

should be taken into consideration.

Your patient hoist should only be used on flat surfaces and with all four casters in contact with the ground.

The patient hoist can be used for transferring in and out of the shower or the bath, but not during showering or bathing.

Riding over obstacles while transferring the patient is completely forbidden.

The patient hoist should not be used as a ladder, nor is it a transport for heavy or hot objects.

When used on mats, carpeted floors or loose floor coverings, the floor covering can get damaged.

Use only Vermeiren approved accessories.


Make sure that the end of the adjustment motor has a safety measurement. Therefore, use the patient hoist only in the middle of the adjustment range and not at the end of the hoist motor.

Compliance with the user and maintenance instructions are an essential part of the guarantee conditions.

Visually impaired people can contact the dealer for the instructions for use.

The patient hoist is suitable for re-use.

1.2 Technical specifications

| Brand | Vermeiren | |
|--|---------------------------------------|---|
| Address | Vermeirenplein 1/15, B-2920 Kalmthout | |
| Type | Patient hoist | |
| Model | Albatros II | |
| | Measurement in drawing | Dimensions |
| Lowest CAP position* | l | 1060 mm |
| Highest CAP position* | k | 1687 mm |
| Lifting range (range of heights) | m | 637 mm |
| Leg length | | 1000 mm |
| Total length | | 1115 mm |
| Total width | | 685 mm |
| Total height | | 1380 mm (actuator lowest position) 1687 mm (actuator highest position) |
| Folded / Dismantled length | | 1115 mm (not foldable) |
| Folded / Dismantled width | | 685 mm (not foldable) |
| Folded / Dismantled height | | 500 mm |
| Minimum leg separation | r | 467 mm |
| Maximum leg separation | q | 769 mm |
| Leg height / Frame height | n | 117 mm |
| Overall width (closed), external dimension | | 584 mm |
| Overall width (opened), castors forwards | P | 945 mm |
| Min. distance between wall / CSP* (minimum height) | w3 (=w2) | 285 mm |
| Min. distance between wall / CSP* (maximum reach) | w2 (=w3) | 285 mm |
| Min. distance between wall / CSP* (maximum height) | w1 | 870 mm |
| Turning circle | | 1160 mm |
| Total weight | | 50,60 kg |
| Weight of chassis + motor + control box | | 34,50 kg |
| Weight of battery | | 2,90 kg |
| Weight of footplate | | 5,10 kg |
| Weight of thorax support | | 1,40 kg |
| Weight of legsupport | | 3,00 kg |
| Weight attachment of legsupport and lifting column | | 3,60 kg |
| Maximum load | | 150 kg |
| Free height, minimum | | 57 mm |
| Maximum reach at 600 mm with legsupport | a1 | 300 mm |
| Maximum reach at 600 mm without legsupport | a2 | 360 mm |
| Maximum reach from support | b | 495 mm |
| Reach from support at separation 700 mm | c | 375 mm |
| Power outlet | | 24V --- max. 250 VA |
| Supply voltage | | 100-240V ~ max. 37-53 VA |
| Maximum current drawn | | max. 400 mA |
| Operating temperature | | +5 to +40°C |
| Relative humidity | | 20% to 90% at 30°C, non-condensing |
| Noise level | | < 50 dB(A) |
| Air pressure | | 700 to 1060 hPa |
| Control box | | Linak CBJ2 |
| Battery | | Linak BAJ1 |
| Manual control | | Linak HB5X (< 5N force required to operate) |
| Motor | | Linak LA34 (7500N) |
| Control box protection class | | IPX4 |

| | | | |
|----------------------------------|---------------------------------------|---|--|
| Brand | Vermeiren | | |
| Address | Vermeirenplein 1/15, B-2920 Kalmthout | | |
| Type | Patient hoist | | |
| Model | Albatros II | | |
| | Measurement in drawing | Dimensions | |
| Battery protection class | | IPX5 | |
| Manual control protection class | | IPX5 | |
| Motor protection class | | IP54 | |
| Insulation class | | II - Type B | |
| Operational performance | | approx. 40 lifts per charge | |
| Periodic operation | | max. 10%, or 2 min. continuous operation/18 min. pauses | |
| Battery capacity | | 2.9 Ah | |
| Emergency stop | | Yes | |
| Manual lowering in emergency | | Yes | |
| Electrical lowering in emergency | | No | |

We reserve the right to introduce technical changes. Measurement tolerance ± 15 mm / 1,5 kg


1 = Highest position, 2 = Maximum reach, 3 = Lowest position

* CSP = Central suspension point

Table 1: Technical specifications

The patient hoist complies to the requirements set up in:
ISO 7176-8: Requirements and test methods for static, impact and fatigue strengths.
ISO 7176-16: Resistance to ignition of upholstered parts.

1.3 Drawing

1.4 Explanation of the symbols

Observe the safety instructions!

Read the instruction manual before use!

Separate recovery and recycling of electric and electronic devices.

CE declaration of conformity

1.5 Delivery

The Vermeiren Albatros II patient hoist shall be delivered with:

- chassis including 4 casters (2 with brakes)
- pole, including handles for pushing
- cantilever arm with attachment for sling
- control box (+ included battery charger)
- battery
- manual control
- motor
- instruction manual
- foot plate + footrests with velcro fastening
- leg support
- thorax support

Before use check if everything is included and that no products are damaged (example by transport, ...).

Please note that the basic configuration may vary from one European country to another. Please contact the specialist dealer in your country for more information.

2 Use

This chapter describes the everyday use. **These instructions are for the user and the specialist dealer.**

The patient hoist is delivered fully assembled by your specialist dealer. The instructions intended for the specialist dealer on how to set up the patient hoist are given in § 3.

2.1 General notes

The patient hoist must only be operated by qualified staff who have been instructed or trained in its specific application.


The use of the patient hoist is restricted to indoors only. Please ensure that it is only used on flat surfaces. It must be used strictly indoors only.

We would like to point out that sources of electromagnetic radiation (e.g. mobile phones, etc.) can cause interference and that the lift's electronics can also affect other electrical appliances.

Even if you have been instructed by your dealer about the operational elements of your lift and their use, we recommend that you read the following pages carefully.

2.1.1 Instructions for lifting patient

Phase I – Preparation. In this phase, the patient should stay on the footplate and put his feet in the special footrests with the velcro fastening. Then stabilize your feet and legs under your knees by using the velcro belt on the foot- and leg support. Make sure that the feet and legs are firmly touching the pillow from the foot- and leg support. Fasten the thorax belt with the four straps to the patient hoist.

Phase II- Lifting. Use the manual control slowly to lift the patient up. The patient must hold the handgrips during lifting.

Phase III – Stabilization. In the last step after lifting to the vertical position, the patient is stabilized by thorax belts. Thorax belts are mounted on both sides of the thorax support. Afterwards you can remove the vest. To complete stabilization you should use the waist belt which stabilizes the hip and trunk of the cross section. In this way, the patient is fully protected in four support points: feet, knees, hips and thorax.

Buttons for waist belt attachment

View properly how the patient is stabilized.

2.2 Use of the JUMBO system (battery, control box + included battery charger)

JUMBO is a modular system combining an actuator, a control box and a battery in a flexible solution, specially developed for patient hoists.

2.2.1 Battery

The battery is above the control box and is replaceable with an integrated clip system.

| | |
|-----------------------------|---------------------------------------|
| Primary voltage | 24 VDC |
| Ambient temperature | +5°C to +40°C |
| Ambient storage temperature | -10° to +50°C |
| Relative storage humidity | max. 90% (non-condensing) |
| Conformity | tested in accordance with IEC 60601-1 |

2.2.2 Control box + included battery charger

The control box is designed with a red emergency stop and an internal battery charger.

| | |
|--------------------------------------|---|
| Primary voltage | 100 - 240 VAC / 50/60 Hz |
| Secondary voltage (charging voltage) | 24 VDC, max. 250 VA |
| Secondary current (charging current) | max. 10 A |
| Safeguards | Protected against reverse polarisation, electrical surges and extreme temperature |
| Ambient temperature | +5°C to +40°C |
| Ambient storage temperature | -10° to +50°C |
| Relative storage humidity | max. 90% (non-condensing) |
| Conformity | tested in accordance with IEC 60601-1 |

We reserve the right to introduce technical changes.

2.2.3 External charger (OPTIONAL)

The external charger shall be mounted on the wall. It is possible to charge the spare battery (optional) with this external charger and you have no waiting time to use the patient hoist when charging the batteries.

| | |
|-----------------------------|---------------------------|
| Primary voltage | 100 - 240 VAC / 50/60 Hz |
| Charging current | max. 650 mA |
| Ambient temperature | +5°C to +40°C |
| Ambient storage temperature | -10° to +50°C |
| Relative storage humidity | max. 90% (non-condensing) |

2.2.4 Charging the batteries

Use only the control box with included battery charger on the Albatros II or the external charger meant to be mounted on the wall.

We recommend that you should regularly recharge the batteries in order to ensure that the hoist remains usable and to extend the lifespan of the batteries. The control box will give a warning sound to tell you when the battery charge is on 50% capacity.

- **FIRST USE**

First, insert the supply cable with the connector side into the appropriate connector on the control box. The plug side of the supply cable must be insert in the socket outlet. The loading time is about 24 hours.

- **RECHARGING**

⚠ WARNING: Risk of injury – Do not use the lift when the supply cable is connected to the socket outlet.

When the batteries are recharging there will burn a yellow light (above on) and an orange light (above charge).

After the recharging process is completed, you should always remove the mains plug from the socket first and then remove the connector plug from the control box.

Do not use the lift while it is being recharged!


If the batteries are not used for a significant period, they discharge slowly by themselves (deep discharge). It then becomes impossible to recharge them with the battery charger supplied. You should therefore recharge the batteries at least once a month even if they are not being used.

2.3 Operating the lift

Please observe the technical information according to which the patient hoists (Albatros II) may be operated. The patient hoist must only be operated by authorised staff who have been trained in its use and operation.

⚠ WARNING: Risk of burns – Be careful when driving in hot or cold environments (sunshine, extreme cold, saunas, etc.) for a sufficient amount of time and when touching - Surfaces can assume environment temperatures.


2.3.1 Altering the chassis width

To get the lift around a wheelchair or other seating furniture, or to increase the stability of the lift in the standing position, the separation of the chassis legs can be increased.


Stand behind the assembled lift and take hold of the handles for pushing the patient hoist (on the left and right, next to the control box). Press the lever at the bottom of the chassis (left or right) down with your foot, and the separation of the legs of the chassis can be increased or decreased.

2.3.2 Parking brakes

Secure the two casters at the back of the chassis by pressing the caster brake plates down gently to their end stops with the tip of your foot. To release the brake, push the brake plate gently back up again with the front of your foot until the casters are free.


2.3.3 Foot plate

The footplate ① can be mounted or removed with the hooks ② to/from the chassis of the patient hoist. If it is required for a greater stabilization of your feet you can use the special footrests ③ with velcro straps for stabilization.

2.3.4 Leg support


⚠ WARNING: Risk of injury – Do not drop the leg support on the patient when adjusting.

The leg support ① can be mounted on the bar of the patient hoist ② with the safety pin ③.

The leg support is meant to keep the legs of the patient in the correct position and to support the legs of the patient during stand up.

The leg supports can be adjusted in depth and 9 different heights.


| Height from footplate to top of leg support | Position |
|---|----------|
| 690 mm | Hole 1 |
| 665 mm | Hole 2 |
| 640 mm | Hole 3 |
| Do not use | Hole 4 |
| Do not use | Hole 5 |
| 495 mm | Hole 6 |
| 470 mm | Hole 7 |
| 445 mm | Hole 8 |
| 420 mm | Hole 9 |

Table 2: Height adjustment leg support

2.3.5 Thorax support


To stabilize and support your thorax, the lift is equipped with thorax pillow and velcro belts to stabilize.

The thorax support can be adjusted in 2 different positions 1045 mm, 1200 mm. These positions are measured from the footplate until the top of the thorax support.

1. Remove the safety pin ①.
2. Adjust the thorax support in the desired position. The position from attachment plate ② shall be like shown in figure.
3. Mount the safety pin ① securely.

2.3.6 Raising and lowering the cantilever arm

The manual control allow the cantilever arm to be adjusted smoothly to any point.

- ① = Raise the cantilever arm
- ② = Lower the cantilever arm

⚠ When using the lift function until the end of the hoist motor the limit switch is activated by electronics (security function). Use the lifting function only in the middle of the adjustment range.

Force required to activate the pushbutton: < 5N

Position the manual control to the top tube of the handles for moving the patient hoist. This is the most ergonomic position.

2.3.7 Moving the patient hoist

⚠ WARNING: Risk of clamping – Be careful passing through restricted passages (e.g. doors).

1. Stand behind the patient hoist and grasp the handles (left and right, next to the battery and control box) with both hands.
2. Make sure that both back caster brakes have been released.
3. Push the patient hoist slowly into the desired position.

When it is not loaded, we advise you to pull the patient hoist backwards, which makes it easier to get round any obstacles (door frames, corners of rooms or furniture).

To move the patient hoist use only the handles for pushing and no other components (actuator, manual control, ...) .

2.4 Slings

⚠ WARNING: Risk of injury – Only use appropriate slings for the patients.

⚠ WARNING: Risk of injury – Do not use damaged slings.

The Albatros II patient hoist must only be used with slings designed for the Vermeiren Albatros II patient hoists.

The applicability and instructions for use for the specific slings should be followed.

2.5 Safety rules

Some safety tips are given below for your security:

- ⚠ Only authorised staff who have been trained in using the lift must use and operate the patient lift.
- ⚠ Make sure when using the lift that there is sufficient room around and above it, as the adjustment movements could otherwise lead to damage or injury.
- ⚠ Always be aware of the lifting arm to avoid injury.
- ⚠ The patient lift must only be used on level surfaces, with all 4 casters contacting the ground evenly.
- ⚠ It must not be used in humid environments.
- ⚠ Only slings that have been designed and approved for the patients must be used (see usage instructions for the various slings). Use of any other slings is at your own risk.
- ⚠ Please follow the instructions of the nursing staff or suitably trained persons, so that no injuries are caused during use of the lift.
- ⚠ Check the medical condition of the patient and be sure that it is possible to lift the patient with this patient hoist. (Refer to intended use)
- ⚠ Be careful with sources of ignition such as lit cigarettes as they may set the sling alight.
- ⚠ Never exceed the maximum load (**150 kg**). The control box is disabled if the lift is overloaded.
- ⚠ When using slings with a maximum load, the lowest load is used between the patients hoist and sling.
- ⚠ The manufacturer accepts no liability for damages and injury caused by inappropriate actions or by not observing the instructions given in this manual.
- ⚠ Use only the handles for pushing to move the patient hoist and no other components (actuator, manual control, ...).
- ⚠ Pull the lifting arm before use manually upwards and check that the lifting arm falls down by himself. If not, contact your specialist dealer.

2.6 Emergency

In cases for emergency the emergency stop and the function for lowering in an emergency must be used.

2.6.1 Emergency stop

In emergencies, and during assembly and disassembly, the emergency stop is activated by pressing the red button on the control box.


Push on the button and turn in the direction of the arrows. If the button is fully pressed to the back turn back in the opposite direction of the arrows.

The emergency stop is deactivated by turning the red button in the direction of the arrows.

2.6.2 Manual emergency lowering

⚠ WARNING: Risk of injury – Adjust the manual emergency lowering for the weight of the patient.

Emergency lowering if the power fails or the batteries are empty can be done using the red pull button (1) at the lower end of the motor. This emergency lowering is factory-set for a patient weight of 50 kg.

Standard calibration: 2500N, 15 mm/s to lower the patient.

Please note that the manual emergency lowering is only possible when the patient is sitting in the patient lift.

Pull force: 10 mm for force from 80 à 100N.

3 Installation and adjustment

The instructions in this chapter are for the specialist dealer.


To find a service facility or specialist dealer near you, contact the nearest Vermeiren facility. A list of Vermeiren facilities can be found on the last page.

⚠ WARNING: Risk of unsafe limitations - Use only the limitations described in this manual.

3.1 Swapping the battery

- We shall decline all liability for damage caused by the use of improper battery.
- Do not use the battery at temperatures below +5°C or above +50°C (the ideal is +20°C).
- If the battery is opened, all liability of the manufacturer shall be voided as well as any claim.

The battery can be swapped by an integrated clip system.


To swap the battery pull on the integrated clip system. The plate will be lifted and the battery can take off the holder (control box).

3.2 Assembly or dismantling


3.2.1 Unpacking

⚠ WARNING: Risk of injury - After unpacking and before further assembly, always first confirm that the emergency stop (the red pushbutton on the control box) is pressed in.

The patient hoist packaging has been selected to provide the best possible protection during transport.


1. Take the patient hoist out of the box and check that all the items have been delivered and that the individual parts do not have any visible defects. If there is any damage, please contact your sales office.


2. Before assembly, put on the caster brakes to make sure that the lift cannot be moved unintentionally. Push the caster brake plate down gently to its end stop with the front of your foot. To release the brake, push the brake plate gently back up again with the front of your foot until the caster is free.

3.2.2 Assembly


⚠ WARNING: Risk of clamping - Take care that no body parts and wires get trapped, crushed or cut during assembly.


1. Remove the safety spindle (locking pin) from the bottom end of the pole by gently pressing the button at the head of the spindle. The safety spindle can now be removed easily.

2. Undo the transport lock (wing screw) gently and the pole can now be raised.


- 2.a** The standing position safety (traction screw) clicks audibly into place once the pole is put in its final position.


3. **⚠ WARNING: Risk of injury - Make sure that the safety spindle is correctly inserted.**


Fix the pole in place by pushing the safety spindle at the end of the pole in with the pushbutton depressed until its other end clearly protrudes. Release the pushbutton on the safety spindle and check that it cannot be pulled out.


4. Remove the cantilever arm from the top end of the pole ②. This is possible by pulling on the pen ①.


5. Mount the slide column ① with the foot support attachment ② in the tube of the patient hoist frame. Retighten the lever ③.


6. Attach the slide and slide holder ① to the cantilever arm ② with the locking pin ③.


7. Pull the actuator gently out the black clip ①.
8. Attach the actuator on the cantilever arm by the locking pin ②. Take care that the actuator is positioned very well between the cantilever arm.

9. Attach the footplate with the hooks to the chassis of the patient hoist.


10. Attach the leg support ① with the locking pin ②. Position the leg support ① to the desired height and depth.


11. Mount the thorax support ① to the slide column ② with the locking pin ③ in a comfortable height for every patient.

3.2.3 Dismantling


- ⚠ **WARNING:** Risk of clamping - Take care that no body parts and wires get trapped, crushed or cut during dismantling.
- ⚠ **WARNING:** Risk of injury - Check the emergency stop (red pushbutton on the control box) is activated, to avoid unintentional adjustment movements.


1. Secure the two casters at the back of the chassis by pressing the caster brake plates down gently to their end stops with the tip of your foot. To release the brake, push the brake plate gently back up again with the front of your foot until the casters are free.


2. Remove the thorax support ① from the slide column ② with the locking pin ③.

3. Remove the leg support ① by the locking pin ②.


4. Remove the footplate from the chassis of the patient hoist.


5. Dismantle the actuator (motor) from the cantilever arm by removing the locking pin ①.


6. Move the actuator to the bottom position and place the actuator in the black clip ②.


7. Remove the slide and slide holder ① from the cantilever arm ② with the locking pin ③.


8. Loosen the lever ①. Remove the slide column ② and foot support attachment ③ out the tube of the patient hoist frame.

9. Move the cantilever arm to the bottom position ①. Fasten the pen on the top end of the pole in the hole of the cantilever arm ②.


10. Remove the safety spindle (locking pin) from the bottom end of the pole by gently pressing the button at the head of the spindle. The safety spindle can now be removed easily.

11. **⚠ WARNING: Risk of injury - Once the pole is unlocked, its own weight can cause considerable downward forces to be exerted.**

Pull out the standing position safety (traction screw) until the pole is released and can be tilted. The pole can now be tilted forwards.

Tip the pole as far as possible forward.

12. **⚠ WARNING: Risk of injury - Make sure that the safety spindle is correctly inserted.**

Fix the pole in place by pushing the safety spindle at the end of the pole in with the pushbutton depressed until its other end clearly protrudes. Release the pushbutton on the safety spindle and check that it cannot be pulled out again.

4 Maintenance

For the maintenance manual of the patient hoists refer to the Vermeiren website: www.vermeiren.be.

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ENGLISH

CONTRACTUAL WARRANTY

We offer 5 years of warranty on standard wheelchairs, lightweight wheelchairs 4 years. Electronic wheelchairs, tricycles, beds and other products: 2 years (batteries 6 months) and multiposition wheelchairs 3 years. This warranty is limited to the replacement of defective or spare parts.

APPLICATION CONDITIONS

In order to claim this warranty, part "B" of this card has to be given to your official Vermeiren dealer. The warranty is only valid when parts are replaced by Vermeiren in Belgium.

EXCEPTIONS

This warranty is not valid in case of:

- damage due to incorrect usage of the wheelchair,
- damage during transport,
- involvement in an accident,
- a dismount, modification or repair carried outside of our company and/or official Vermeiren dealership,
- normal wear of the wheelchair,
- non-return of the warranty card

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All rights reserved, including translation.
Vermeiren Suisse S.A.
Vermeiren Nederland B.V.