## BM 57

<b>(D)</b>	<b>Blutdruckmessgerät</b> Gebrauchsanweisung.....	2-17
<b>(GB)</b>	<b>Blood pressure monitor</b> Instructions for use.....	18-32
<b>(F)</b>	<b>Tensiomètre</b> Mode d'emploi .....	33-48
<b>(E)</b>	<b>Tensiómetro</b> Instrucciones de uso .....	49-64
<b>(I)</b>	<b>Misuratore di pressione</b> Istruzioni per l'uso .....	65-79
<b>(TR)</b>	<b>Tansiyon ölçme cihazı</b> Kullanım Kılavuzu .....	80-94
<b>(RUS)</b>	<b>Прибор для измерения кровяного давления в плечевой артерии</b> Инструкция по применению .....	95-112
<b>(PL)</b>	<b>Ciśnieniomierz</b> Instrukcja obsługi .....	113-129
	<b>Electromagnetic Compatibility Information.....</b>	<b>130</b>

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## Dear customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage, beauty and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

With kind regards,  
Your Beurer team

## 1. Getting to know your device

Check that the packaging of the Beurer BM 57 blood pressure monitor has not been tampered with and make sure that all the required contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

The upper arm blood pressure monitor is used to carry out non-invasive measurement and monitoring of the arterial blood pressure values in adults.

This allows you to quickly and easily measure your blood pressure, save the measured values and display the development and average values of the measured values taken.

You are also warned of possible existing cardiac arrhythmia. The recorded values are classified and evaluated graphically.

Store these instructions for use for future reference and make them accessible to other users.



## 2. Important notes

### Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and the accessories:

	Attention
	Note Note on important information
	Observe the instructions for use
	Application part, type BF
	Direct current
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Manufacturer

 <p>Storage -20°C 55°C RH ≤85%</p>	Permissible storage temperature and humidity
 <p>Operating 10°C 40°C RH ≤85%</p>	Permissible operating temperature and humidity
IP 22	Protected against ingress of solid foreign objects greater than 12.5 mm in diameter. Protected against drops of water falling at up to 15° from vertical.
SN	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

### Notes on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measured values taken by you are for your information only – they are no substitute for a medical examination. Discuss the measured values with your doctor and never base

any medical decisions on them (e.g. medicines and their administration).

- Do not use the blood pressure monitor on newborns or patients with preeclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the device on people who have the specified upper arm measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions remove the cuff from the arm.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or intravascular therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.

- Do not place the cuff over wounds as this may cause further injury.
- Please note that data transfer and data storage is only possible when your blood pressure monitor is supplied with power. As soon as the batteries are empty, the blood pressure monitor loses the date and time.
- To conserve the batteries, the blood pressure monitor switches off automatically if you do not press any buttons for 30 seconds.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

### **Instructions for storage and maintenance**

- The blood pressure monitor is made from precision and electronic components. The accuracy of the measured values and service life of the device depend on its careful handling:
  - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
  - Do not drop the device.
  - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
  - Only use the cuff included with the delivery or original replacement parts. Otherwise incorrect measured values will be recorded.
- Do not press the buttons before the cuff is placed on the arm.

### **Notes on handling batteries**

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- Choking hazard! Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- Risk of explosion! Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush the batteries.

### **Battery disposal**

- The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
- The codes below are printed on batteries containing harmful substances:  
 Pb = Battery contains lead,  
 Cd = Battery contains cadmium,  
 Hg = Battery contains mercury


### Instructions for repairs and disposal

- Batteries do not belong in household waste. Please dispose of empty batteries at the collection points intended for this purpose.
- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorised suppliers.
- Do not open the device. Failure to comply will invalidate the warranty.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the device at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.


### 3. Device description

1. Cuff
2. Cuff line
3. Cuff connector
4. Display
5. Memory buttons **M1/M2**
6. START/STOP button 
7. Scale for classifying the measurements
8. Connection for cuff connector (left-hand side)
9. Connection for mains part

## Information on the display:

1. Time/date
2. Systolic pressure
3. Diastolic pressure
4. Calculated pulse value
5. Pulse symbol ♥
6. Release air (arrow)
7. Number of memory space/memory display for average value ( $\bar{P}$ ), morning ( $P_M$ ), evening ( $P_E$ )
8. Classification of measurements
9. Battery display symbol
10. Cardiac arrhythmia symbol ♥
11. User memory  $P_1$   $P_2$
12. Symbol for *Bluetooth*® transfer

## System requirements for the Beurer “HealthManager” PC software.

### 1. Supported operating systems:

- Windows Vista SP2
- Windows 7
- Windows 7 SP1
- Windows 8
- Windows 10

### 2. Supported architectures:

- x86 (32 bit)
- x64 (64 bit)


### 3. Hardware requirements:

- Recommended: at least Pentium 1 GHz or faster with at least 1 GB RAM
- Free memory on the primary partition of at least:
  - x86 – 600 MB
  - x64 – 1.5 GB
- Graphic resolution from: 1024 x 768 pixels

## 4. Preparing the measurement

### Inserting the batteries

- Remove the battery compartment lid on the rear of the device.
- Insert 2x 1.5 V AA (alkaline type LR06) batteries. Make sure that the batteries are inserted the correct way round. Do not use rechargeable batteries.
- Close the battery compartment lid again carefully.

All display elements are briefly displayed, 24 h flashes on the display. Now set the date and time as described below.

If the battery replacement symbol  is permanently displayed, you can no longer perform any measurements and must replace all batteries. Once the batteries have been removed from the device, the date and time must be set again. Any saved measured values are retained.

### Set the hour format, date, time and Bluetooth® settings

The following section describes the functions and settings available on the blood pressure monitor.

**Hour format** → **Date** → **Time** → **Bluetooth®**

It is essential to set the correct date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

- If you press and hold the **M1** or **M2** memory button, you can set the values more quickly.

Press and hold the **START/STOP**  button for 5 seconds if you have not reinserted the batteries.

#### Hour format

The hour format now flashes on the display.

- Select the desired hour format using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The year flashes on the display.

- Select the year using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The month flashes on the display.

- Select the month using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The day flashes on the display.

- Select the current day using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

- If you have set the 12h hour format, the month is displayed before the day.

The hours flash on the display.

- Select the current hours using the **M1/M2** memory buttons and confirm with the **START/STOP** button **(I)**.

The minutes flash on the display.

- Select the current minutes using the **M1/M2** memory buttons and confirm with the **START/STOP** button **(I)**.

The *Bluetooth*® symbol flashes on the display.

- Use the **M1/M2** memory buttons to select whether automatic *Bluetooth*® data transfer is activated (*Bluetooth*® symbol flashes) or deactivated (*Bluetooth*® symbol is not shown) and confirm with the **START/STOP** button **(I)**.

**(i)** *Bluetooth*® transfers will reduce the battery life.

## Operation with the mains part

You can also operate this device with a mains part. When doing so, there must not be any batteries in the battery compartment. The mains part can be obtained from specialist retailers or from the service address using order number 071.51.

- To prevent possible damage to the device, the blood pressure monitor must only be used with the mains part described here.
- Insert the mains part into the connection provided for this purpose on the right-hand side of the blood pressure monitor. The mains part must only be connected to the mains voltage that is specified on the type plate.

- Then insert the mains plug of the mains part into the mains socket.
- After using the blood pressure monitor, unplug the mains part from the mains socket first and then disconnect it from the blood pressure monitor. As soon as you unplug the mains part, the blood pressure monitor loses the date and time setting but the saved measured values are retained.

## 5. Measuring blood pressure


Ensure the device is at room temperature before measuring. The measurement can be performed on the left or right arm.

### Attaching the cuff

Place the cuff on to the bare left upper arm. The circulation of the arm must not be hindered by tight clothing or similar.

The cuff must be placed on the upper arm so that the bottom edge is positioned 2 – 3 cm above the elbow and over the artery. The line should point to the centre of the palm.

Now tighten the free end of the cuff, but make sure that it is not too tight around the arm and close the hook-and-loop fastener. The cuff should be fastened so that two fingers fit under the cuff.

Now insert the cuff line into the connection for the cuff connector.


- ❗ If the measurement is performed on the right upper arm, the line should be located on the inside of your elbow. Ensure that your arm is not pressing on the line.

Blood pressure may vary between the right and left arm, which may mean that the measured blood pressure values are different. Always perform the measurement on the same arm.

If the values between the two arms are significantly different, please consult your doctor to determine which arm should be used for the measurement.

**Important:** The unit may only be operated with the original cuff. The cuff is suitable for an arm circumference of 23 to 43 cm.

### Adopting the correct posture

- Before every measurement, relax for about five minutes. Otherwise deviations can occur.
- You can take the measurement while sitting or lying down. Always make sure that the cuff is at heart level.
- To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on

something. Do not cross your legs. Place your feet flat on the ground.

- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.

### Performing the blood pressure measurement

As described above, attach the cuff and adopt the posture in which you want to perform the measurement.

- Press the **START/STOP** button ❶ to start the blood pressure monitor. All display elements are briefly displayed.


The blood pressure monitor will begin the measurement automatically after 3 seconds. The measurement is taken during the inflation phase.

- ❗ You can cancel the measurement at any time by pressing the **START/STOP** button ❶.

As soon as a pulse is found, the pulse symbol ♥ will be displayed.

- The systolic pressure, diastolic pressure and pulse rate measurements are displayed.

- EE appears if the measurement could not be performed properly (see chapter 10 "Error messages/troubleshooting"). Repeat the measurement.

- Now select the desired user memory by pressing the **M1** or **M2** memory buttons. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant  $r_1$  or  $r_2$  symbol appears on the display.
- Press the **START/STOP** button **I** to switch off the blood pressure monitor. The measurement is then stored in the selected user memory.
- If the device is not switched off manually, it will switch off automatically after approx. 3 minutes.

If **Bluetooth® data transfer** has been activated, data is transferred after having confirmed the user memory and pressing the **START/STOP** button **I**.

- The **Bluetooth®** symbol flashes on the display. The blood pressure monitor now attempts to connect to the app for approx. 30 seconds.
- The **Bluetooth®** symbol stops flashing as soon as a connection is established. All measurement data is automatically transferred to the app. Once the data has been successfully transferred, the device switches off. If the data transfer was unsuccessful, “EE 7” appears on the display.
- If a connection to the app cannot be established after 30 seconds, the **Bluetooth®** symbol goes out and the blood pressure monitor switches off automatically after 3 minutes.

- ⓘ Please note that you must add the blood pressure monitor in “My devices” in the Beurer “HealthManager” app to enable data transfers. The Beurer “HealthManager” app must be active to allow data transfers.  
If the latest data is not displayed on your smartphone, repeat the data transfer as described in chapter 8.

If you forget to turn off the blood pressure monitor, it will switch off automatically after a further 30 seconds. In this case too, the value is stored in the selected or most recent user memory and the data is transferred if **Bluetooth®** data transfer has been activated.

- Wait at least 5 minutes before taking another measurement.

## 6. Evaluating results

### Cardiac arrhythmia:

This device can identify potential disruptions of the heart rhythm when measuring and if necessary, indicates this after the measurement with the symbol .

This can be an indicator for arrhythmia. Arrhythmia is a condition in which the heart rhythm is abnormal because of flaws in the bioelectrical system that regulates the heartbeat. The symptoms (skipped or premature heart beats, pulse being slow or too fast) can be caused by factors such as heart disease, age, physical make-up, excess stimulants, stress or lack of sleep.

Arrhythmia can only be determined through an examination by your doctor.

If the symbol  is shown on the display after the measurement has been taken, repeat the measurement. Please ensure that you rest for 5 minutes beforehand and do not speak or move during the measurement. If the symbol  appears frequently, please consult your doctor.

Self-diagnosis and treatment based on the measurements can be dangerous. Always follow your GP's instructions.

### **Classification of measurements:**

The measurements can be classified and evaluated in accordance with the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The classification on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the values of systole and diastole fall into two different categories (e.g. systole in the 'High normal' category and diastole in the 'Normal' category), the graphical classification on the device always shows the higher category; for the example given this would be 'High normal'.

<b>Blood pressure value category</b>	<b>Systole (in mmHg)</b>	<b>Diastole (in mmHg)</b>	<b>Action</b>
Setting 3: severe hypertension	≥ 180	≥ 110	seek medical attention
Setting 2: moderate hypertension	160–179	100–109	seek medical attention
Setting 1: mild hypertension	140–159	90–99	regular monitoring by doctor
High normal	130–139	85–89	regular monitoring by doctor
Normal	120–129	80–84	self-monitoring
Optimal	< 120	< 80	self-monitoring

Source: WHO, 1999 (World Health Organization)

## 7. Display and delete measured values

The results of every successful measurement are stored together with the date and time. The oldest measurement is overwritten in the event of more than 60 measurements.

- Select the desired user memory (M1 M2) using the **M1** or **M2** memory buttons if the device is switched off.

– To view the measurement data for user memory M1, press the **M1** memory button.

– To view the measurement data for user memory M2, press the **M2** memory button.

The average of all measurements appears on the display.

If *Bluetooth*® is activated (the symbol will flash on the display), the blood pressure monitor is attempting to connect to the app. If you press the **M1** button in the meantime, the transfer will be cancelled and the morning average will be displayed.

If you press the **M2** button, the transfer is cancelled and the measurement data of user memory M2 is displayed. The  symbol is no longer shown. The buttons are deactivated as soon as a connection is established and the data is transferred.

- ① Press the **M1** button if you have selected user memory 1.  
Press the **M2** memory button if you have selected user memory 2.

**M** flashes on the display.

The average value of all saved measured values in this user memory is displayed.

- Press the relevant memory button (**M1** or **M2**).

**M** flashes on the display.

The average value of the morning measurements for the last 7 days is displayed (morning: 5.00 a.m. – 9.00 a.m.).

Press the relevant memory button (**M1** or **M2**).

**M** flashes on the display.

The average value of the evening measurements for the last 7 days is displayed (evening: 6.00 p.m. – 8.00 p.m.).

When the relevant memory button (**M1** or **M2**) is pressed again, the last individual measurement is displayed (in this example, measurement 03).

When the relevant memory button (**M1** or **M2**) is pressed again, you can view your individual measurements.

- To switch the device off again, press the START/STOP button **ⓘ**.

**ⓘ** You can exit the menu at any time by pressing the **START/STOP** button **ⓘ**.

- To clear the memory of the relevant user memory, you must first select a user memory.
- Start the retrieval of the average measured values. **R** flashes on the display and the average value of all saved measured values in this user memory is displayed.
- Then press and hold the **M1** and **M2** memory buttons for 5 seconds.

All the values in the current user memory are deleted.

To clear individual measurements from the relevant user memory, you must first select a user memory.

- Start the retrieval of the individual measured values.
- Press and hold the memory buttons **M1** and **M2** for 5 seconds.
- The selected value is deleted. The device briefly displays **CL 00**.
- If you would like to delete other values, repeat the process described above.

You can switch off the device at any time by pressing the START/STOP button.

### Transfer via **Bluetooth® Smart**

It is also possible to transfer the measured values saved on the device to your smartphone using **Bluetooth® Smart**.

You will need the Beurer “HealthManager” app for this. The app is available from the App Store.

Proceed as follows to transfer values:

If **Bluetooth®** is activated in the settings menu, the data is transferred automatically after the measurement. The **ⓘ** symbol appears in the top left of the display (see chapter 4 “Preparing the measurement”).

### Step 1: BM 57

Activate *Bluetooth*® on your device (see chapter 4 “Preparing the measurement, *Bluetooth*®”).

### Step 2: Beurer “HealthManager” app

In the Beurer “HealthManager” app, add the BM 57 under “Settings/My devices”.

### Step 3: BM 57

take a measurement.

### Step 4 BM 57:

data transfer immediately following measurements.

- If *Bluetooth*® data transfer has been activated, data is transferred after having confirmed the user memory and pressing the **START/STOP** button .

### Step 4 BM 57:

Data transfer at a later point:

- Go to memory mode (chapter 7). Select the desired user memory. The *Bluetooth*® transfer starts automatically.

 The Beurer “HealthManager” app must be active to allow data transfers.

If your smartphone has a protective cover, remove this to ensure that there is no interference during the transfer. Begin the data transfer in the Beurer “HealthManager” app.

## 8. Cleaning and storing the unit

- Clean the device and cuff carefully using a slightly damp cloth only.
- Do not use any cleaning agents or solvents.
- Under no circumstances should you hold the unit under water, as this can cause liquid to enter and damage the unit.
- If you store the device, no heavy objects should be placed on top of it. The cuff line should not be bent sharply.

In such cases, repeat the measurement and/or data transfer. Ensure that the cuff line is correctly attached and that you do not move or speak during the measurement.

The *Bluetooth*® word mark and logos are registered trademarks owned by *Bluetooth*® SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

## 9. Error messages/troubleshooting

In the event of errors, the error message **EE** appears on the display.

Error messages may appear if:

- Systolic or diastolic pressure could not be measured (**EE 1**)
- The cuff is faulty or fastened too tightly or too loosely (**EE 2**)
- The pump pressure is higher than 290 mmHg (**EE B**)
- The data could not be sent via *Bluetooth*® (**EE 7**)

## 10. Technical specifications

Model No.	BM 57
Measurement method	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement range	Cuff pressure 0 – 300 mmHg, Systolic 50 – 250 mmHg, Diastolic 40 – 180 mmHg, Pulse 40– 180 beats/minute
Display accuracy	Systolic $\pm 3$ mmHg, Diastolic $\pm 3$ mmHg, Pulse $\pm 5$ % of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: Systolic 8 mmHg / Diastolic 8 mmHg
Memory	2 x 60 memory spaces
Dimensions	L 94 mm x W 161 mm x H 64 mm
Weight	317 g (without cuff)
Cuff size	23 to 43 cm
Permissible operating conditions	10°C to +40°C, $\leq 85\%$ relative humidity (non-condensing)
Permissible storage conditions	-20 to +55°C, $\leq 85\%$ relative humidity, 800– 1050 hPa ambient pressure
Power supply	2x 1.5 V  AA batteries

Battery life	Lasts for approx. 150 measurements depending on the blood pressure and pump pressure as well as the number of <i>Bluetooth</i> <sup>®</sup> connections.
Accessories	Cuff, instructions for use, storage pouch
Classification	Internal supply, IP22, no AP or APG, continuous operation, application part type BF

Data transfer via Bluetooth<sup>®</sup> wireless technology

The blood pressure monitor uses *Bluetooth*<sup>®</sup> Smart (low energy), 2.4 GHz frequency band, Compatible with *Bluetooth* 4.0 smartphones/tablets

List of supported smartphones/tablets

Technical information is subject to change without notification to allow for updates.

- This device complies with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Custom-

er Services address or found at the end of the instructions for use.

- This device corresponds to the EU Medical Devices Directive 93/42/EEC, the German Medical Devices Act (Medizinproduktegesetz) and the standards DIN EN ISO81060-1 (non-invasive sphygmomanometers – Part 1: Requirements and test methods for non-automated measurement type). EN1060-3 (non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC80601-2-30 (Medical electrical equipment – Part 2 – 30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.
- We hereby guarantee that this product complies with the European R&TTE Directive 1999/5/EC. Please contact the specified service address to obtain more detailed information such as the CE conformity declaration.

## 11. Mains part

Model no.	8194047/02
Input	100–240 V AC, 50–60 Hz
Output	6 V, 600 mA, in conjunction with Beurer blood pressure monitors only
Manufacturer	Ansmann AG
Protection	The device is double protected and has a primary-side cutout switch which disconnects the device from the mains in case of malfunction. Ensure that you have removed the batteries from the battery compartment before you use the mains part.
	Insulated/protection class 2
Housing and protective covers	The housing of the mains part protects users from touching live parts or parts that could be live (for example with finger, needle, checking hook). Do not touch the patient and the output connector of the AC/DC mains part at the same time.

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Table 1

### For all ME EQUIPMENT and ME SYSTEMS

#### *Guidance and manufacture's declaration – electromagnetic emission*

The BM57 is intended for use in the electromagnetic environment specified below. The customer of the user of the BM57 should assure that it is used in such and environment.

<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The BM57 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The BM57 Blood Pressure Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

**Table 2**  
**For all ME EQUIPMENT and ME SYSTEMS**

*Guidance and manufacture's declaration – electromagnetic immunity*

The BM57 is intended for use in the electromagnetic environment specified below. The customer or the user of BM57 should assure that it is used in such an environment.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. The BM57 can continue the operation during power mains interruptions due to the usage of battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Mains power quality should be that of a typical commercial or hospital environment.
<b>NOTE:</b> $U_T$ is the a.c. mains voltage prior to application of the test level.			

**Table 3**  
**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

*Guidance and manufacture's declaration – electromagnetic immunity*

The BM57 is intended for use in the electromagnetic environment specified below. The customer or the user of BM57 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the BM57, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b></p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BM57 is used exceeds the applicable RF compliance level above, the BM57 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the BM57.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Table 4**  
**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

*Recommended separation distances between portable and mobile RF communications equipment and the BM 85*

The BM57 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BM57 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BM57 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.