SANITAS

SBM 37

D	Blutdruckmessgerät	
	<i>Gebrauchsanweisung</i>	2-16
GB	Blood pressure monitor	
	<i>Instructions for use</i>	17-31
F	Tensiomètre	
	<i>Mode d'emploi</i>	32-46
I	Misuratore di pressione	
	<i>Istruzioni per l'uso</i>	47-61
RUS	Прибор для измерения кровяного давления в плечевой артерии	
	<i>Инструкция по применению</i>	62-78
	Electromagnetic Compatibility Information	79-82

CE 0483

Serviceadresse:
Hans Dinslage GmbH
Riedlinger Straße 28
88524 Uttenweiler, GERMANY
Tel.-Nr.: +49 (0) 73 74 / 91 57 66
E-mail: service@sanitas-online.de

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ENGLISH

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

1. Getting to know your device

Check that the packaging of the Sanitas SBM 37 blood pressure monitor has not been tampered with and make sure that the required contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

The upper arm blood pressure monitor is used to carry out non-invasive measurement and monitoring of the arterial blood pressure values in adults.

This allows you to quickly and easily measure your blood pressure, save the measured values and display the development and average values of the measured values taken.


You are also warned of possible existing cardiac arrhythmia. The recorded values are classified and evaluated graphically. Store these instructions for use for future reference and make them accessible to other users.

2. Important notes

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and the accessories:

	Attention
	Note Note on important information
	Observe the instructions for use
	Application part, type BF
	Direct current
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Manufacturer

	Permissible storage temperature and humidity
	Permissible operating temperature and humidity
	Protect from moisture
SN	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

Notes on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measured values taken by you are for your information only – they are no substitute for a medical examination. Discuss the measured values with your doctor and never base

any medical decisions on them (e.g. medicines and their administration).

- Do not use the blood pressure monitor on newborns or patients with preeclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the device on people who have the specified upper arm measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions remove the cuff from the arm.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or intravascular therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.

- Do not place the cuff over wounds as this may cause further injury.
- Please note that data transfer and data storage is only possible when your blood pressure monitor is supplied with power. As soon as the batteries are empty, the blood pressure monitor loses the date and time.
- To conserve the batteries, the blood pressure monitor switches off automatically if no buttons are pressed for 30 seconds.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Instructions for storage and maintenance

- The blood pressure monitor is made from precision and electronic components. The accuracy of the measured values and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.
 - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
 - Only use the cuff included with the delivery or original replacement parts. Otherwise incorrect measured values will be recorded.
- Do not press the buttons before the cuff is placed on the arm.

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- Choking hazard! Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- Risk of explosion! Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush the batteries.

Battery disposal

- The empty, completely discharged batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
- The codes below are printed on batteries containing harmful substances:
 - Pb = Battery contains lead,
 - Cd = Battery contains cadmium,
 - Hg = Battery contains mercury


i Instructions for repairs and disposal

- Batteries do not belong in household waste. Please dispose of empty batteries at the collection points intended for this purpose.
- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorised suppliers.
- Do not open the device. Failure to comply will invalidate the warranty.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the device at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.


3. Device description

1. Cuff
2. Cuff line
3. Cuff connector
4. Display
5. Memory buttons **M1/M2**
6. START/STOP button **I**
7. Scale for classifying the measurements
8. Connection for cuff connector (left-hand side)

Information on the display:

1. Time/date
2. Systolic pressure
3. Diastolic pressure
4. Calculated pulse value
5. Pulse symbol
6. Release air (arrow)
7. Number of memory space/memory display for average value (M), morning (M_M), evening (M_E)
8. Classification of measurements
9. Battery display symbol
10. Cardiac arrhythmia symbol
11. User memory
12. Symbol for *Bluetooth*[®] transfer

System requirements for the HealthCoach app

Bluetooth[®]4.0, iOS from Version 7.0, Android[™] devices from Version 4.3 with Bluetooth[®] Smart Ready


List of compatible devices:

4. Preparing the measurement

Inserting the batteries

- Remove the battery compartment lid on the rear of the device.
- Insert four 1.5 V AAA (alkaline type LR03) batteries. Make sure that the batteries are inserted the correct way round. Do not use rechargeable batteries.
- Close the battery compartment lid again carefully.

All display elements are briefly displayed, 24 h or 12 h flashes in the display. Now set the date and time as described below.

If the battery replacement symbol is permanently displayed, you can no longer perform any measurements and must replace all batteries. Once the batteries have been removed from the device, the date and time must be set again. Any saved measured values are retained.

Set the hour format, date, time and **Bluetooth**[®] settings

The following section describes the functions and settings available on the blood pressure monitor.

Hour format

Date

Time

Bluetooth[®]

It is essential to set the correct date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

i If you press and hold the **M1** or **M2** memory button, you can set the values more quickly.

Press and hold the START/STOP button ① for 5 seconds.

Hour format

The hour format now flashes on the display.

- Select the desired hour format using the **M1/M2** memory buttons and confirm with the **START/STOP** button ①.

Date

The year flashes on the display.

- Select the year using the **M1/M2** memory buttons and confirm with the **START/STOP** button ①.

The month flashes on the display.

- Select the month using the **M1/M2** memory buttons and confirm with the **START/STOP** button ①.

The day flashes on the display.

- Select the current day using the **M1/M2** memory buttons and confirm with the **START/STOP** button ①.

- i** If you have set the 12h hour format, the month is displayed before the day.

Time

The hours flash on the display.

- Select the current hours using the **M1/M2** memory buttons and confirm with the **START/STOP** button ①.

The minutes flash on the display.

- Select the current minutes using the **M1/M2** memory buttons and confirm with the **START/STOP** button ①.

The *Bluetooth*® symbol flashes on the display.

- Use the **M1/M2** memory buttons to select whether automatic *Bluetooth*® data transfer is activated (*Bluetooth*® symbol flashes) or deactivated (*Bluetooth*® symbol is not shown) and confirm with the **START/STOP** button **ⓘ**.

ⓘ *Bluetooth*® transfers will reduce the battery life.

5. Measuring blood pressure




Ensure the device is at room temperature before measuring. The measurement can be performed on the left or right arm.

Attaching the cuff


Place the cuff on to the bare left upper arm. The circulation of the arm must not be hindered by tight clothing or similar.

The cuff must be placed on the upper arm so that the bottom edge is positioned 2 – 3 cm above the elbow and over the artery. The line should point to the centre of the palm.

Now tighten the free end of the cuff, but make sure that it is not too tight around the arm and close the hook-and-loop fastener. The cuff should be fastened so that two fingers fit under the cuff.

Now insert the cuff line into the connection for the cuff connector.

ⓘ If the measurement is performed on the right upper arm, the line should be located on the inside of your elbow. Ensure that your arm is not pressing on the line.


Blood pressure may vary between the right and left arm, which may mean that the measured blood pressure values are different. Always perform the measurement on the same arm.

If the values between the two arms are significantly different, please consult your doctor to determine which arm should be used for the measurement.

Important: The unit may only be operated with the original cuff. The cuff is suitable for an arm circumference of 22 to 36 cm.

A larger cuff for upper-arm circumferences of 35 to 44 cm can be obtained from specialist retailers or from the service address using order number 163.387.

Adopting the correct posture

- Before every measurement, relax for about five minutes. Otherwise deviations can occur.
- You can take the measurement while sitting or lying down. Always make sure that the cuff is at heart level.
- To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground.
- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.

Performing the blood pressure measurement

As described above, attach the cuff and adopt the posture in which you want to perform the measurement.

- Press the **START/STOP** button **(i)** to start the blood pressure monitor.
All display elements are briefly displayed.


The blood pressure monitor will begin the measurement automatically after 3 seconds.

- You can cancel the measurement at any time by pressing the **START/STOP** button **(i)**.

As soon as a pulse is found, the pulse symbol **♥** will be displayed.

- The systolic pressure, diastolic pressure and pulse rate measurements are displayed.

- **Er** appears if the measurement could not be performed properly (see chapter 9 "Error messages/troubleshooting"). Repeat the measurement.

Er

- Now select the desired user memory by pressing the **M1** or **M2** memory buttons. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant r_1 or r_2 symbol appears on the display.
- Press the **START/STOP** button **ⓘ** to switch off the blood pressure monitor. The measurement is then stored in the selected user memory.
- If the device is not switched off manually, it will switch off automatically after 3 minutes.

If **Bluetooth® data transfer** has been activated, data is by transferred after having confirmed the user memory by pressing the pressing the **START/STOP** button **ⓘ**.

- The **Bluetooth®** symbol on the display flashes and the blue LED lights up. The blood pressure monitor now attempts to connect to the app for approx. 30 seconds.
- The **Bluetooth®** symbol stops flashing as soon as a connection is established. All measurement data is automatically transferred to the app. Once the data has been successfully transferred, the device switches off. If the data transfer was unsuccessful, the blue LED goes out and "Er 7" appears on the display.
- If a connection to the app cannot be established after 30 seconds, the **Bluetooth®** symbol goes out and the blood pressure monitor switches off automatically after 3 minutes.

- ⓘ Please note that you must add the blood pressure monitor in "My devices" in the "HealthCoach" app to enable data transfers. The "HealthCoach" app must be active to allow data transfer. If the latest data is not displayed on your smart-phone, repeat the data transfer as described in chapter 7.

If you forget to turn off the blood pressure monitor, it will switch off automatically after approximately 3 minutes. In this case too, the value is stored in the selected or most recent user memory and the data is transferred if **Bluetooth®** data transfer has been activated.

- Wait at least 5 minutes before taking another measurement.

6. Evaluating results

Cardiac arrhythmia:

This device can identify potential disruptions of the heart rhythm when measuring and if necessary, indicates this after the measurement with the symbol .

This can be an indicator for arrhythmia. Arrhythmia is a condition in which the heart rhythm is abnormal because of flaws in the bioelectrical system that regulates the heartbeat. The symptoms (skipped or premature heart beats, pulse being slow or too fast) can be caused by factors such as heart disease, age, physical make-up, excess stimulants, stress or lack of sleep.

Arrhythmia can only be determined through an examination by your doctor.

If the symbol  is shown on the display after the measurement has been taken, repeat the measurement. Please ensure that you rest for 5 minutes beforehand and do not speak or move during the measurement. If the symbol  appears frequently, please consult your doctor.

Self-diagnosis and treatment based on the measurements can be dangerous. Always follow your GP's instructions.

Classification of measurements:

The measurements can be classified and evaluated in accordance with the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The classification on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the values of systole and diastole fall into two different categories (e.g. systole in the 'High normal' category and diastole in the 'Normal' category), the graphical classification on the device always shows the higher category; for the example given this would be 'High normal'.

Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
Setting 3: severe hypertension	≥ 180	≥ 110	seek medical attention
Setting 2: moderate hypertension	160–179	100–109	seek medical attention
Setting 1: mild hypertension	140–159	90–99	regular monitoring by doctor
High normal	130–139	85–89	regular monitoring by doctor
Normal	120–129	80–84	self-monitoring
Optimal	< 120	< 80	self-monitoring

Source: WHO, 1999 (World Health Organization)

7. Display and delete measured values

The results of every successful measurement are stored together with the date and time. The oldest measurement is overwritten in the event of more than 60 measurements.

- Select the desired user memory (M1 M2) with the **M1** or **M2** memory buttons if the device is switched off.

– To view the measurement data for user memory M1, press the **M1** memory button.

– To view the measurement data for user memory M2, press the **M2** memory button. The average of all measurements appears on the display.

If *Bluetooth*® is activated (the  symbol flashes on the display), the blood pressure monitor attempts to connect to the app. If you press the **M1** button in the meantime, the transfer is cancelled and the average values are displayed. If you press the **M2** button, the transfer is cancelled and the measurement data of user memory M2 is displayed. The  symbol is no longer shown. The buttons are deactivated as soon as a connection is established and the data is transferred.

- ① Press the **M1** button if you have selected user memory 1.

Press the **M2** memory button if you have selected user memory 2.

A flashes on the display.

The average value of all saved measured values in this user memory is displayed.

- Press the relevant memory button (**M1** or **M2**).

M flashes on the display.

The average value of the morning measurements for the last 7 days is displayed (morning: 5.00 a.m. – 9.00 a.m.).

Press the relevant memory button (**M1** or **M2**).

P flashes on the display.

The average value of the evening measurements for the last 7 days is displayed (evening: 6.00 p.m. – 8.00 p.m.).

When the relevant memory button (**M1** or **M2**) is pressed again, the last individual measurement is displayed (in this example, measurement 03).

When the relevant memory button (**M1** or **M2**) is pressed again, you can view your individual measurements.

- To switch the device off again, press the **START/STOP** button **ⓘ**.

ⓘ You can exit the menu at any time by pressing the **START/STOP** button **ⓘ**.

- To clear the memory of the relevant user memory, you must first select a user memory.
- Start the retrieval of the average measured values. **Ⓜ** flashes on the display and the average value of all saved measured values in this user memory is displayed.
- Press and hold the memory button **M1** or **M2** for 5 seconds depending on the user memory you are in.

All the values in the current user memory are deleted.

To clear individual measurements from the relevant user memory, you must first select a user memory.

- Start the retrieval of the individual measured values.
- Press and hold the memory button **M1** or **M2** for 5 seconds (depending on the user memory you are in).
- The selected value is deleted. The device briefly displays **CL 00**.
- If you would like to delete other values, repeat the process described above.

You can switch off the device at any time by pressing the **START/STOP** **ⓘ** button.

Transfer via **Bluetooth® Smart**

It is also possible to transfer the measured values saved on the device to your smartphone using **Bluetooth® Smart**.

You will need the HealthCoach app for this. The app is available from the App Store.

Proceed as follows to transfer values:

If **Bluetooth®** is activated in the settings menu, the data is transferred automatically after the measurement. The **Ⓜ** symbol appears in the top left of the display (see chapter 4 „Preparing the measurement“).

Step 1: SBM 37

Activate *Bluetooth*® on your device (see chapter 4 "Preparing the measurement", *Bluetooth*®).

Step 2: "HealthCoach" app

In the HealthCoach app, add the SBM 37 under "Settings/My devices".

Step 3: SBM 37

take a measurement.

Step 4: SBM 37

data transfer immediately following measurements.

- If *Bluetooth*® data transfer has been activated, data is by transferred after having confirmed the user memory by pressing the **START/STOP** button ①.

Step 4: SBM 37

Data transfer at a later point:

- Go to memory mode (chapter 7). Select the desired user memory. The *Bluetooth*® transfer starts automatically.

① The HealthCoach app must be active to allow data transfer. If your smartphone has a protective cover, remove this to ensure that there is no interference during the transfer. Begin the data transfer in the HealthCoach app.

8. Cleaning and storing the unit

- Clean the device and cuff carefully using a slightly damp cloth only.
- Do not use any cleaning agents or solvents.
- Under no circumstances should you hold the unit under water, as this can cause liquid to enter and damage the unit.
- If you store the device, no heavy objects should be placed on top of it. The cuff line should not be bent sharply.

In such cases, repeat the measurement and/or data transfer. Ensure that the cuff line is correctly attached and that you do not move or speak during the measurement.

9. Error messages/troubleshooting

In the event of errors, the error message E_r appears on the display.

Error messages may appear if:

- it was not possible to record the pulse: E_r1 ;
- you move or speak during the measurement: E_r2 ;
- the cuff is fastened too tightly or loosely: E_r3 ;
- errors occur during the measurement: E_r4 ;
- the pump pressure is higher than 300 mmHg: E_r5 ;
- the batteries are almost empty : E_r6 ;
- E_r7 the data could not be sent via *Bluetooth*®.

The word mark *Bluetooth*[®] and accompanying logo are registered trademarks of *Bluetooth*[®]SIG, Inc. Any and all use of these marks by Hans Dinslage GmbH is done so under licence. Other trademarks and trade names are the property of the relevant holder.

10. Technical specifications

Model No.	SBM 37
Measurement method	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement range	Cuff pressure 0 – 250 mmHg, Systolic 50 – 250 mmHg, Diastolic 30 – 200 mmHg, Pulse 40– 180 beats/minute
Display accuracy	Systolic ± 3 mmHg, Diastolic ± 3 mmHg, Pulse ± 5 % of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: Systolic 8 mmHg / Diastolic 8 mmHg
Memory	2 x 60 memory spaces
Dimensions	L 139 mm x W 94 mm x H 48 mm
Weight	236 (without cuff)
Cuff size	22 to 36 cm
Permissible operating conditions	+5°C to +40°C, 15-93% relative humidity (non-condensing)

Permissible storage conditions	-25°C to +70°C, $\leq 93\%$ relative humidity, 800– 1050 hPa ambient pressure
Power supply	4x 1.5 V \equiv AAA batteries
Battery life	Lasts for approx. 200 measurements depending on the blood pressure and pump pressure as well as the number of <i>Bluetooth</i> [®] connections.
Accessories	Cuff, instructions for use, storage pouch
Classification	Internal supply, IPX0, no AP or APG, continuous operation, application part type BF
Data transfer via Bluetooth [®] wireless technology	The blood pressure monitor uses <i>Bluetooth</i> [®] Smart (low energy), 2.4 GHz frequency band, compatible with <i>Bluetooth</i> 4.0 smartphones/tablets
	List of supported smartphones/tablets

Technical information is subject to change without notification to allow for updates.

- This device complies with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and

mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Services address or found at the end of the instructions for use.

- The device complies with the EU Medical Devices Directive 93/42/EC, the German Medical Devices Act (Medizinproduktegesetz) and the standards EN1060-1 (non-invasive sphygmomanometers – Part 1: General requirements), EN1060-3 (non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC80601-2-30 (Medical electrical equipment – Part 2 – 30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.
- We hereby guarantee that this product complies with the European R&TTE Directive 1999/5/EC. Please contact the specified service address to obtain more detailed information such as the CE conformity declaration.

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Electromagnetic compatibility

Electromagnetic emissions

Guidance and manufacturer's declaration – electromagnetic emissions for all EQUIPMENT and SYSTEMS

The Blood pressure monitor (SBM 37) is intended for use in the electromagnetic environment specified below. The customer of the user of the Blood pressure monitor (SBM 37) should assure that it is used in such and environment.

Emission test	Compliance	Electromagnetic environment–guidance
RF emissions CISPR 11	Group 1	The Blood pressure monitor (SBM 37) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The Blood pressure monitor (SBM 37) is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Electromagnetic immunity

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

The Blood pressure monitor (SBM 37) is intended for use in the electromagnetic environment specified below. The customer of the user of the Blood pressure monitor (SBM 37) should assure that it is used in such and environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment–guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.

Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
--	------	------	---

Guidance and manufacturer's declaration – electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

The Blood pressure monitor (SBM 37) is intended for use in the electromagnetic environment specified below. The customer or the user of Blood pressure monitor (SBM 37) should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Blood pressure monitor (SBM 37), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance.</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.5 GHz</p>

			<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.</p> <p>To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Blood pressure monitor (SBM 37) is used exceeds the applicable RF compliance level above, the Blood pressure monitor (SBM 37) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Blood pressure monitor (SBM 37).</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for EQUIPMENT or SYSTEM that are not LIFESUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the Blood pressure monitor (SBM 37).

The Blood pressure monitor (SBM 37) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Blood pressure monitor (SBM 37) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Blood pressure monitor (SBM 37) as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0.1167	0.2334
0,1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3786
100	11.6667	23.3334

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.