# HAND EXERCISES

IN THE FOLLOWING LANGUAGES:

**ENGLISH • DUTCH • FRENCH • GERMAN • ITALIAN • SPANISH**

Daily use conditions and strengthens the hands, wrists and forearms by flexing and toning the muscles, ligaments and tendons. These exercises are designed to do with any of our hand exercises.

Begin all exercises slowly and use low resistance & few repetitions. Hold each position 3-5 seconds and relax - repeat 5 to 10 times. Little by little increase to 3 sets of 10 (30 reps). When graduating to a higher resistance begin again with 5 to 10 reps and build slowly.

If you experience pain or fatigue – stop immediately. As with any exercise, excessive or incorrect use can lead to pain or injury. Therefore progress slowly and increase reps and resistance very cautiously.

## EN- Hook Position

Place finger tips on individual buttons and position ergonomic palm bar with the hook over the web space between the thumb and index finger. Keeping finger tips slightly bent, flex all fingers toward the center of the unit.

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## EN- Tip-to- Tip Pinch

Place the tip of the thumb on the center of the ergonomic palm bar. Position the finger to be exercised on one of the center buttons. Keeping both tips slightly bent, flex the thumb and the finger being exercised toward the center of the unit at the same time.

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## EN- Gross Grasp

Place the middle of the fingers on the buttons with the finger tips wrapped over the top. Positioning the hook of the ergonomic palm bar over the web space between the thumb and index finger, flex the thumb and all fingers toward the center of the unit as if making a fist.

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## EN- Fingertip Flexion

Place all finger tips on their own button, with the tip of the thumb on the center of the ergonomic base. Keeping all finger tips bent, flex the thumb and all the fingers toward the center of the unit at the same time.

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## **EN- Key Pinch**

Place the hook of the ergonomic palm bar on the inside of the bent index finger. With the tip of the thumb slightly bent, flex the thumb on the first button.

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## **EN- Trigger Pinch**

Loosely hold the unit vertically with the buttons toward the fingers and the hook of the ergonomic palm bar positioned over the web space between the thumb and index finger. With the middle of the index finger, press the second button as if squeezing a trigger.

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