

beurer  
**CONNECT**



UNIVERSAL  
Cuff



Download on the  
**App Store**

ANDROID APP ON  
**Google play**




Data Privacy/  
Data Security \*

www.tuv.com  
ID 000042517



HealthManager



Kostenloser Zugriff auf Beurer HealthManager Free access to Beurer HealthManager

medical

beurer

## BM 57 Tensiomètre Bluetooth®

Mesure entièrement automatique de la pression sanguine et du pouls au niveau du bras

Connexion innovante entre le Smartphone et le tensiomètre. Mesures suivies via Beurer HealthManager App / Web / PC

Données transférées via Bluetooth® Smart Technologie

Mesure automatique par inflation: pour une mesure rapide et confortable

Rétroéclairage blanc LCD

2 x 60 emplacements de mémoire/ Moyenne des valeurs

Classification des mesures sur une échelle de valeurs

Détection de l'arythmie: met en garde en cas de turbulence du rythme cardiaque

Produit médical

Contrôle optimal au niveau du bras pour les tours de bras entre 23 - 43 cm

Date et heure / Arrêt automatique

Message en cas d'erreur d'utilisation

Indicateur de changement des piles

Passeport de pression artérielle

Adaptateur disponible séparément Art-Nr.: 071.51

2 piles AAA de 1,5V incluses

3 ans de garantie

UE: 4 / Carton de transport: 2 x 4

EAN-Nr.: 4211125 65822 9

Art.-Nr.: 658.22

## BM 57 Bluetooth® Blood pressure monitor

Fully automatic blood pressure and pulse measurement on the upper arm

Innovative networking between smartphone and blood pressure monitor. Measurements can be monitored using Beurer HealthManager app/Internet/PC

Transfer using Bluetooth® Smart technology

Automatic measurement during inflation: for fast and comfortable measurement

LCD with white backlighting

2 x 60 memory spaces / Average value calculation

Classification of measurement with color chart

Arrhythmia detection: warning in case of possible heart rhythm disturbance

Medical product

Optimal control of the upper arm circumference with cuff for upper arm circumferences 23 - 43 cm

Date and time / automatic switch off

Signal in case of application error

Low battery indicator

Blood pressure passport

Main part available separately  
Item no.: 071.51

Inclusive 2 x 1.5V batteries AAA

3 year guarantee

Sales unit: 4 / Shipping carton: 2 x 4

EAN no.: 4211125 65822 9

Item no.: 658.22

\*Application en ligne vérifiée Beurer HealthManager cloud

\* audited online application Beurer HealthManager cloud