



Dispo  
Octobre 2016  
PVC: 49,99€



\_ avant le repas

\_ après le repas

\_ marquage général



**Codefree**

## GL 43 mg/dL Lecteur de glycémie connecté

Facile à utiliser : grand écran, bandes de test larges, marquage facile de la valeur mesurée pour une prise en charge optimisée du diabète : avant le repas, après le repas, marquage général

Dépistage de sous-alimentation : mesure uniquement sur un échantillon sanguin suffisamment grand

Kit de démarrage : appareil de massage, autopiqueur, 5 lancettes, 5 bandes de test, piles, étui, mode d'emploi, manuel abrégé, journal de glycémie sanguine, câble USB

Adapté à l'auto-évaluation

0,6 µL de sang suffisant

Temps de mesure : env. 5 secondes

Sans code

480 valeurs mémorisables avec date et heure (sauvegardées même lors d'un changement de piles)

Les résultats sauvegardés peuvent défilés de haut en bas

Valeurs moyennes sur 7, 14, 30 ou 90 jours pour toutes les mesures et pour les valeurs marquées

Échantillon de sang : sang complet frais (capillaire)

Etalonnage plasma

Plage de mesure: 20-630 mg/dL

Avec câble USB et téléchargement gratuit du logiciel HealthManager sur: [www.beurer.com](http://www.beurer.com) > SERVICE > Téléchargement

Abschaltautomatik

Arrêt automatique

1 pile lithium CR 2032 pour plus de 1 000 mesures

Dimensions (LxH): 88 x 51 x 15 mm

Poids: env. 41 g (piles incluses)

Stockage/Transport:	Temp.:	Humidité
Appareil:	-10°C à +65°C	< 90 %
Bandes de test:	+2°C à +30°C	< 90 %
En fonction:	+10°C à +40°C	< 90 %

Garantie 3 ans

UE: 5 / Carton de transport: 4x5

EAN-Nr.: 4211125 46318 2

Art.-Nr.: 463.18



## GL 43 mg/dL Blood Glucose Monitoring System

Easy to use: Extra large display, wide test strips, simple measurement marking for optimised diabetes management: before eating, after eating, general marking

Blood volume check: measurement only possible with sufficient sample of blood

Starter Pack: measuring unit, lancing device, 5 lancets, 5 test strips, batteries, pouch, instruction manual, short instruction, blood glucose diary, USB cable

Suitable for self-testing

Only 0.6 µl of blood necessary

Period of measurement: approx. 5 seconds

Codefree

480 stored results with date, time (data memory also with battery change)

Stored results can be scrolled forward and back

Average value calculation for 7, 14, 30 or 90 days for - all measurements - marked values

Blood sample: fresh whole blood (capillary)

Plasma calibrated

Measurement range: 20-630 mg/dL

With USB cable and HealthManager software download for free under: [www.beurer.com](http://www.beurer.com) > SERVICE > download

Automatic switch off

Medical product

1 x Lithium battery CR 2032 for more than 1.000 measurements

Product measurements (LxWxH): 88 x 51 x 15 mm

Product weight: approx. 41 g (including batteries)

Storage/Transport:	Temp.:	Rel.humidity:
Unit:	-10°C to +65°C	< 90 %
Test stripes:	+2°C to +30°C	< 90 %
Operating:	+10°C to +40°C	< 90 %

3 year guarantee

Sales unit: 5 / Shipping carton: 4x5

EAN no.: 4211125 46318 2

Item no.: 463.18

