




nicht final not final



Data Privacy/
Data Security *
www.tuv.com
ID 000042517



Kostenloser Zugriff auf Beurer HealthManager Free access to Beurer HealthManager

german|engineering

medical

PO 60 Bluetooth® Oxymètre de pouls

Mesure la saturation pulsée en oxygène (SpO₂) et la fréquence cardiaque (pouls)

Connexion innovante entre le Smartphone et l'oxymètre de pouls. Suivi des mesures via Beurer HealthManager App / Web / PC

Transmission des données par Bluetooth® Smart 4.0

100 espaces de mémoire inclus

Mesure simple et sans douleur

Petit et léger - pour une utilisation partout

Particulièrement adapté aux personnes souffrant de:

- Insuffisance cardiaque
- Brancho-pneumopathie chronique obstructive
- Asthme bronchique

Peut également être utilisé dans le cadre de sports de haute altitude (par ex. alpinisme, ski, aviation sportive)

Ecran couleur très lisible

Ecran avec 6 affichages disponibles

Luminosité de l'écran réglable

Affichage graphique du pouls

Indicateur du niveau des piles

Arrêt automatique

Inclus: bandoulière, pochette pour la ceinture

Produit médical

2 piles AAA incluses

Dimensions (L/W/H): 52 x 30 x 33 mm

Poids: env. 27 g (sans les piles)

Garantie 3 ans

UE: 4 / Carton de transport: 8 x 4

EAN-Nr.: 4211125 45420 3

Art.-Nr.: 454.20

*Application en ligne certifiée BeurerHealthManager Cloud

PO 60 Bluetooth® Pulse Oximeter

Determining arterial oxygen saturation (SpO₂) and heart rate (pulse frequency)

Innovative networking between smartphone and pulse oximeter. Measurements can be monitored using Beurer HealthManager app/Internet/PC

Transfer using Bluetooth® Smart technology

Incl. 100 memory spaces

Simple and completely pain-free measurement

Small and light – for use at home and on the move

Particularly suitable for persons with:

- Heart failure
- Chronic obstructive pulmonary diseases
- Bronchial asthma

Suitable for sports at high altitudes (e.g. mountaineering, skiing and aviation sports)

Easy-to-read colour display

Display with 6 available views

Adjustable display brightness

Graphical display of pulse

Low battery indicator

Automatic switch-off

Incl. retaining strap and belt pouch

Medical product

Incl. 2 x AAA batteries

Product measurements (L/W/H): 52 x 30 x 33 mm

Product weight: approx. 27 g (without batteries)

3 year guarantee

Sales unit: 4 / Shipping carton: 2 x 4

EAN no.: 4211125 45420 3

Item no.: 454.20

* audited online application Beurer HealthManager cloud

beurer

0615 Irrtum und Änderungen vorbehalten