



## The short VACOPed For Foot and Ankle Injuries

VACOPed Short is the shorter and lighter version of VACOPed. Thanks to the integrated VACO12 technology, the orthosis adapts perfectly to the foot and ankle. VACOPed Short allows not only immobilisation of the ankle, but also a defined range of motion (ROM) from 15° of dorsal extension to 20° of plantar flexion.

### Your Advantages at a Glance

- Excellent pressure relief and high stability through VACO12 technology
- Optimum adaptation to individual anatomical requirements
- Functional after-treatment possible due to defined ROM from -15° to +20° (adjustable in 5° increments)
- Easy wound control simplifies treatment
- Improved hygiene through changeable and washable liners

### Indications

- Ankle joint fractures
- Metatarsal fractures
- Achilles tendon ruptures
- Calcaneal fractures
- Severe sprains
- Ligament rupture
- Ligament surgery
- Arthrodesis of the metatarsal joint
- Hammer and claw toes
- Fracture of the hallux



- 1 Vacuum cushion **VACO12** individual adjustment, even with swelling
- 2 Pull-out toe section for maximum protection of the forefoot
- 3 Strong, high frame provides stability
- 4 Adjustable joint  
ROM: -15° to +20° (adjustable in 5° increments)  
**Immobilisation**  
Allows fixation of the ankle joint  
**Mobilisation**  
Use as a functional orthosis
- 5 Changeable sole/Step-in-sole  
Rocker sole for a physiological gait
- 6 Partial load control: **Clickdisk**

### Sizes

Available in 3 sizes, universal for left or right foot.

Size	Small	Medium	Large
EUR	35 - 40	40 - 45	45 - 50
UK	2 ½ - 6 ½	6 ½ - 10 ½	10 ½ - 14
US (M)	3 ½ - 7 ½	7 ½ - 11 ½	11 ½ - 15
US (W)	4 ½ - 8 ½	8 ½ - 12 ½	12 ½ - 16